

Product Information

ZmTech® Improved quantitative PCR Mixture

Catalog Number: **A207010**

Storage Temperature: **-20°C**

Kit contains:

An Improved qPCR Mixture (2x), containing 2x PCR reaction Buffer, Heat-Activated Taq DNA Polymerase, dNTPs (dATP, dCTP, dGTP, dTTP), 6mM MgCl₂, PCR Buffer Stabilizers, RNase/DNase I inhibitors, Green I dye chemistry and Internal reference.

Description:

ZmTech® Improved quantitative PCR Mixture is a 2x reaction mix designed for the quantitative real-time PCR assays where template is limited or very low-copy-number targets. It has been developed to minimize primer-dimers, which are of particular concern when using dye chemistry for detection.

ZmTech® Improved qPCR Mixture contains low levels of internal reference dye, ideal for those block-based platforms which are more sensitive to dye levels. This internal reference dye is necessary for ABI sequence detection systems, and no interfere with LightCycler® I/II systems.

Suggestions for use: The Taq DNA polymerase is inactive at room temperature and a 10-minute hot start at 95°C is essential to active the enzyme. Hot start procedure can significantly improve the results of DNA amplifications by reducing the generation of nonspecific amplification products and primer-dimer artifacts.

Storage and stability:

Store this kit at -20°C in a constant-temperature freezer. Avoid prolonged exposure to light.

Reaction conditions/suggestions for use:

- RNA preparation: Integrity and Purity of RNA template is essential. Impurities (proteins, phenol, chloroform) may interfere the qPCR Reactions and detection systems.
- MgCl₂ Optimization: MgCl₂ is a co-factor for the Taq DNA Polymerase. A final MgCl₂ concentration of 3mM should be ideal for most applications. However, the MgCl₂ requirements for the polymerase often vary, depending on the particular template and primers used. MgCl₂ can be optimized within 3mM-7mM range.
- cDNA template concentration: An excess of template will inhibit PCR reactions. 100ng of cDNA template in a 25ul reaction is normally sufficient. However, this may vary, depending on the particular gene of interest.

Precautions and Disclaimer:

This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

Procedure:

I. Preparation of PCR Master Mix for a single reaction (total volume: 25uL) in a 0.2 ml or 0.5ml microtube

| Component | Volume (μL) | Final Concentration |
|-------------------------------|-------------|----------------------|
| 1. Improved qPCR Mixture (2X) | 12.5 | (1x) |
| 2. MgCl ₂ | * | (3mM-7mM) |
| 3. DNA Template | 2-5 | (determined by user) |
| 4. Forward Primer (5μM) | 1 | 200nM |
| 5. Reverse Primer (5μM) | 1 | 200nM |
| 6. PCR grade water | * | up to 25 μL |

II. Setup typical thermal cycling parameters

| | | |
|----------------------------|--|----------------|
| 1. Enzyme activation Step: | 95°C | 10 minutes |
| 2. 25-40 Cycles: | | |
| Denaturation | 95°C | 15 seconds |
| Annealing | X°C | 30 seconds |
| | Dependent on T _m of primers | |
| Extension | 72°C | (acquire data) |
| | 60 base pairs/second | |

III. Troubleshooting

| | Possible Cause | Recommendation |
|-----------------------------|--|--|
| 1. Low/no Signal | • Too few cycles • Too high annealing temperature • Poor template • Primers degraded | Increase cycles Decrease temperature in increments of 2-4°C Increase template con. Redesign/optimization |
| 2. Multiple Products | • Too high template • Too high primer • Too low annealing temperature • Primer-dimers | Decrease template con. Decrease primer con. in 2 fold dilutions Increase temperature in increments of 2-3°C Hot-start 95°C 10mins |
| 3. No linear standard curve | • Too high or low template • Primer-dimers | Decrease or increase template concentration Hot-start 95°C 10mins |
| 4. Fluo- rescence variable | • Too high or low template • qPCR instrument contaminated • No well mix | Decrease or increase template concentration Decontaminate qPCR instrument Mix qPCR mixture well |