



® • **ZmTech Universal qPCR Mixtures (Cat. K209090)**

**Kit contains:**

Stored at -20°C

- 2x Ready-to-use universal qPCR mixture, including optimization of 2x PCR Buffer, Heat-Activated Taq DNA Polymerase, dNTPs(dATP, dCTP, dGTP, dTTP), 6mM MgCl<sub>2</sub>, PCR Buffer Stabilizer, RNase/DNase I Inhibitors and Internal Reference.

**Reagents required but not provided:**

- PCR grade water and Reverse Transcriptase (RT).

**Description:**

- ZmTech universal qPCR mixture (K209090) is a 2x reaction mixture suitable for most real-time PCR assays, designed for qPCR assays where template is limited or very low-copy-number targets. It has been developed to minimize primer-dimers, which are of particular concern when using fluorescent probe for detection.
- ZmTech universal qPCR mixture contains low levels of internal reference, ideal for those block-based platforms which are more sensitive to dye levels. This internal reference dye is necessary for ABI Sequence detection systems, and no interfere with LightCycler I/II systems.
- **Suggestions for use:** The Taq DNA polymerase is inactive at room temperature and a 5-minute hot start at 95°C is essential to active the enzyme. Hot start procedure can significantly improve the results of DNA amplifications by reducing the generation of nonspecific amplification products and primer-dimer artifacts.

**Reaction conditions/suggestions for use:**

- **RNA preparation:** Integrity and Purity of RNA template is essential. Impurities (proteins, phenol, chloroform) may interfere the RT-PCR Reactions and detection systems.
- **MgCl<sub>2</sub> Optimization:** MgCl<sub>2</sub> is a co-factor for the Taq DNA Polymerase. A final MgCl<sub>2</sub> concentration of 3mM should be ideal for most applications. However, the MgCl<sub>2</sub> requirements for the polymerase often vary, depending on the particular template and primers used. MgCl<sub>2</sub> can be optimized within 3mM-7mM range.
- **Primer Concentration:** A final primer concentration of 200nM is sufficient for most reactions. However, primer concentration can be optimized within 25-900nM range.
- **cDNA Template Concentration:** 50ng of human cDNA in a 25ul reaction is ideal for real-time PCR reactions. However, this may vary, depending on the particular gene of interest. Too high template concentration will inhibit PCR reactions.
- **Probe Concentration:** In most case, using a combination of primers and probe at a molar ratio of 3:1 is sufficient. However, probe concentration can be optimized within the 10-250nM range.



## Procedure: (suggested protocol)

### 1. Preparation of PCR Master Mix for a single reaction (total volume: 25uL) in a PCR microtube or plate

Component	Volume (μL)	Final Concentration
All-IN-One qPCR Mixture(2x)	12.5	1x
MgCl <sub>2</sub>	X	(3mM-7mM)
DNA Template	2	(100pg-500pg/μL)
Forward primer (5μM)	1	200nM
Reverse primer (5μM)	1	200nM
Probe	1	70nM
PCR grade water	up to 25 μL	

### 2. Setup typical thermal cycling parameters

Enzyme activation Step:	95°C	5 minutes	
25-40 Cycles:			
Denaturation	95°C	15 seconds	
Annealing	X°C	30 seconds dependent on Tm of primers	
Extension	60°C	30 seconds (acquire data at the end of this step)	
Or: <u>Run reactions in a 25ul glass capillary tube.</u>			
Step:	Temperature:	Time:	Temperature Transition Rate
1. Denaturation (1 cycle)	94°C	5 minutes	20°C/second
2. Cycles (35 cycles):			
Denature	94°C	15 seconds (none)	20°C/second
Annealing	X°C	15 seconds (none)	20°C/second
Extension	72°C	8 seconds (none)	20°C/second
Detection	72°C	0 seconds (single)	20°C/second
3. Melting(1 cycle)	X°C	1 minute (none)	20°C/second
		(95°C, 0 Second, 0.15°C/second, con.)	
4. Cooling (1cycle)	40°C	30 seconds (none)	20°C/second

### Precautions and Disclaimer:

- This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.