Technical Bulletin ZmTech Scientifique: Innovation and Development of Laboratory Accessories and Bio-Technologies

**Product Information:** 

2x SybrGreen qPCR Master-mix (Cat.: Q2100A/ Q2100B/ Q2100N)

Storage:@-20°C

The 2x qPCR Master Mix (<u>Q2100A</u>) is a ready-to-use hot-start 2x qPCR cocktail containing everything (except a DNA template and primer set), specifically optimized for efficient amplification of GC-rich and/or problematic templates by qPCR. Suitable for qPCR instruments: <u>ABI-7000/7300/7700/7900</u>, <u>Step-One/+</u>.

The 2x qPCR Master Mix (<u>Q2100B</u>) is a ready-to-use hot-start 2x qPCR cocktail containing everything (except a DNA template and primer set), specifically optimized for efficient amplification of GC-rich and/or problematic templates by qPCR. Suitable for qPCR instruments: <u>ABI-7500/Q6/Q7Flex/Vii7</u>, <u>Stratagene MX3000/3005P</u>, <u>Qiagen Rotor-gene 3000</u>.

The 2x qPCR Master Mix (<u>Q2100N</u>) is a ready-to-use hot-start 2x qPCR cocktail containing everything (except a DNA template and primer set), specifically optimized for efficient amplification of GC-rich and/or problematic templates by qPCR. Suitable for qPCR instruments: <u>Roche LightCycler 480/96</u>, <u>MJ-C4/O2 Takara TP-8000</u>, <u>Bio-Rad iQ/iQ5/CFX96/C1000</u>, <u>Qiagen Rotor-gene 6000/G/Q</u>, <u>Pikoreal 96</u>.

## **Suggested PCR Protocol:**

I. Preparation of qPCR Master Mix for a single reaction (total volume: 20uL):

Component	Volume (µL)	Final Concentration	
2x qPCR Mastermix	10	1x 250nM 250nM (1pg-0.5ug) gDNA	
Forward primer (10µM)	0.5		
Reverse primer (10μM)	0.5		
DNA Template	1-5		
PCR grade water	up to 20 μL		

## II. Setup typical thermal cycling parameters

2. 37-45 cycles:			
Denatura	tion 94°C	5 seconds	
Annealin	g 55-60°C	35 seconds (Fluorescent signal	collection)
3. Dissociation stage			

**Precautions and Disclaimer:** This product and procedure described are intended for R&D use only. Purchase of this product does not convey a license to perform any patented process.

ZmTech Scientific endeavours to assist clients based on the highest level of customer service, competitive pricing and customer satisfaction. Our mission is: Convenience, Speed, Safety and Economy.